



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tape umbilical,3mmx50m,non-ster

Product Code :

LBNY-0021-2200067



www.laboratoryinstrumentindia.com

Description :

Non sterile umbilical tape roll.
Material: 100 % cotton.
Ligature thread for umbilical cord.
Braided and flat.
Non-absorbable.
Size selected: Tape umbilical roll.
Length: 50m.
Width: 3mm.
Non-sterile.
Single use.

Technical Specification :

Weight/Volume
Estimated volume: 0.387cdm
Estimated weight: 0.068kg



Laboratory Instrument India

