



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tape umbilical,3mmx50m,non-ster

**Product Code :**

LBNY-0021-2200067



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Non sterile umbilical tape roll.  
Material: 100 % cotton.  
Ligature thread for umbilical cord.  
Braided and flat.  
Non-absorbable.  
Size selected: Tape umbilical roll.  
Length: 50m.  
Width: 3mm.  
Non-sterile.  
Single use.

**Technical Specification :**

Weight/Volume  
Estimated volume: 0.387cdm  
Estimated weight: 0.068kg



**Laboratory Instrument India**

