



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Compress,paraffin,10x10cm,ster/BOX-10

**Product Code :**

LBNY-0021-2200068



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Non-detectable by X-ray.

Material:

Netting material: absorbent gauze, 100% cotton.

Paraffin substance: mixture of balsam of Peru (approximately 1g) and soft paraffin Q. Suff. (approximately 100g).

Components:

Elastic-netting material with large mesh, impregnated with soft paraffin based material.

Sterile.

Disposable.

Initial sterilisation method:

Ethylene oxide gas or ionising radiation.

Size selected: Compress, gauze, paraffin.

Width: approximately 10cm.

Length: approximately 10cm.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.468cdm

Estimated weight: 0.140kg

---

# Laboratory Instrument India

