



Email : sales@laboratoryinstrumentindia.com

Product Name :

Compress, gauze, 10x10cm, ster/PAC-5

Product Code :

LBNY-0021-2200069



www.laboratoryinstrumentindia.com

Description :

Gauze compress, surgical folding, i.e. so that there are no free threads apparent after folding or when the first outside fold is opened. Non-detectable by X-ray.

Thread count warp 95 to 105 threads/dm;

Material: Absorbent gauze, 100% cotton.

Sterile.

Disposable.

Type of gauze: 17 threads/cm² (grammage 23g/m²).

Components: Bleached, purified textile, plain weave.

Number of folds (thickness): 12.

Initial sterilisation method: Ethylene oxide gas or ionising radiation.

Size selected: Compress, gauze, folded.

Width: approximately 10cm.

Length: approximately 10cm

Weft 66 to 74 threads/dm

Weight: Approximately 23g/m².

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.040cdm

Estimated weight: 0.020kg



Laboratory Instrument India