



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrial Transformer Trainer

Product Code :
LBNY-0002-620004



Description :

Electrial Transformer Trainer

Technical Specification :

Transformer provided will be high quality copper wound
Primary will be 440V and secondary of the transformer will be 220V.
Panel consist of MS powder coated housing with back side removable.
Transformer will be differently housed in a MS powder cabinet for easy carrying and other general use.
Electrical panel fitted with front panel of high class insulated Bakelite sheet 6mm thick.
It is fitted with appropriate analog or digital voltmeter, ammeter to record the readings.
Short circuit protection with the help of MCB.

Experiments:-

Single Phase Paralleling Transformer
Single Phase Open Circuit of Transformer
Three Phase Transformer in Parallel
Three Phase Banking of 3, Single Phase
Transformer Technical Features;
Variable Three Phase Power Supply AC
Power supply provides variable Three Phase Voltage

Electrial Transformer Trainer, Electrial Transformer Trainer Equipments, Electrial Transformer Trainer Tools, Electrial Transformer Trainer Tool Kits, Electrial Transformer Trainer Manufacturers, Electrial Transformer

Trainer Suppliers from India, China, Kenya



Laboratory Instrument India