



Email : sales@laboratoryinstrumentindia.com

Product Name :

Infusion giving set,sterile,s.u.

Product Code :

LBNY-0021-220071



www.laboratoryinstrumentindia.com

Description :

Sterile single use device used for parenteral administration of infusions (IV fluids) by gravity.
Precision flow regulator: Smooth roller clamp to facilitate easy, safe control and adjustment of fluid rates.
Drop counting chamber: Transparent drip chamber, calibrated at 20 drops/ml, with 15-20µm fluid filter.
Tubing: Transparent tubing, minimum length 150cm, with injection site latex or latex-free (or Y-injection port), with distal connector preferably Luer Lock connector.
Single use
Sterile
Device components:
Perforator: Sharp piercing device with protective cap.
Air inlet: With integrated bacteriological filter.

Technical Specification :

Infusion giving set, sterile, single use



Laboratory Instrument India

