



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Lancet,blood,ster,disp/PAC-200

**Product Code :**

LBNY-0021-2200072



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sharp point, all sides tapered provides quick, clean puncture ribbed surface for firm grip

Sturdy, ready to use manual blood lancet

Material: Chrome-polished stainless steel

Typical incision depth, approximately 2mm

Single use

Sterile individually packed

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

Dimensions:

Length, approximately 40mm

Width, approximately 6mm

Thickness, approximately 0.05mm.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.336cdm

Estimated weight: 0.100kg.

---

# Laboratory Instrument India

