



Email : sales@laboratoryinstrumentindia.com

Product Name :

Lancet,blood,ster,disp/PAC-200

Product Code :

LBNY-0021-2200072



www.laboratoryinstrumentindia.com

Description :

Sharp point, all sides tapered provides quick, clean puncture ribbed surface for firm grip

Sturdy, ready to use manual blood lancet

Material: Chrome-polished stainless steel

Typical incision depth, approximately 2mm

Single use

Sterile individually packed

Initial sterilisation method: Ethylene oxide gas or Gamma radiation.

Dimensions:

Length, approximately 40mm

Width, approximately 6mm

Thickness, approximately 0.05mm.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.336cdm

Estimated weight: 0.100kg.

Laboratory Instrument India

