



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sut,abs,DEC2,need 3/8 26mm,tri/BOX-36

Product Code :

LBNY-0021-2200073



www.laboratoryinstrumentindia.com

Description :

Type of thread: Absorbable, synthetic braided thread (PGA, polyglycolic acid).

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC 2 (3/0).

Length: 75cm.

Sterile.

Disposable.

Type of needle:

Length: 26mm.

Point: triangular.

Shape: curved needle, 3/8 circle.

Initial sterilisation method:

Ethylene oxide gas.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 1.176cdm

Estimated weight: 0.200kg

Laboratory Instrument India

