



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Sut,abs,DEC3,need 1/2 30mm,round/BOX-36

**Product Code :**

LBNY-0021-2200075



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Type of thread: Absorbable, synthetic braided thread (PGA, polyglycolic acid).

Length: 75cm.

Suture, absorbable, synthetic, braided DEC3 (2/0), curved needle, ½ circle, 30mm, round point, box of 36

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC3 (2/0).

**Technical Specification :**

Sterile.

Type of needle:

Shape: curved needle, 1/2 circle.

Length: 30mm.

Point: round.

Disposable.

Initial sterilisation method:

Ethylene oxide gas.

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg



**Laboratory Instrument India**