



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sut,abs,DEC3,spool/BOX-36

Product Code :

LBNY-0021-2200076



www.laboratoryinstrumentindia.com

Description :

Initial sterilisation method:

Ethylene oxide gas.

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC3 (2/0).

Thread without needle.

Type of thread: Absorbable, synthetic braided thread (PGA, polyglycolic acid).

Length: 150cm.

Sterile.

Disposable.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg



Laboratory Instrument India

