



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Sut,abs,DEC4,need 3/8 50mm,round/BOX-36

**Product Code :**

LBNY-0021-220077



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC4 (1).

Type of thread: Absorbable, synthetic braided thread (PGA, polyglycolic acid).

Length: 75cm.

Type of needle:

Length: Approximately 50mm.

Shape: curved needle, 3/8 circle.

Point: round.

Sterile.

Disposable.

Initial sterilisation method:

Ethylene oxide gas.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg

---

# Laboratory Instrument India

