



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sut,abs,DEC4,need 1/2 30mm,round/BOX-36

Product Code :

LBNY-0021-2200078



www.laboratoryinstrumentindia.com

Description :

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC4 (1).

Type of thread: Absorbable, synthetic braided thread (PGA, polyglycolic acid).

Length: 75cm.

Type of needle:

Point: round.

Shape: curved needle, 1/2 circle.

Length: 30mm.

Sterile.

Disposable.

Initial sterilisation method:

Ethylene oxide gas.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg

Laboratory Instrument India

