



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sut,nonabs,DEC3,need 3/8 30mm,tri/BOX-36

Product Code :

LBNY-0021-2200079



www.laboratoryinstrumentindia.com

Description :

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC3 (2/0).

Type of thread: Non absorbable, synthetic monofilament (Nylon).

Length: 75cm.

Type of needle:

Point: triangular, conventional cutting

Shape: curved needle, 3/8 circle.

Length: approximately 30mm.

Sterile.

Disposable.

Initial sterilisation method:

Ethylene oxide gas.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.030 cdm

Estimated weight: 0.010 kg.

Laboratory Instrument India

