



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sut,nonabs,DEC1,need 1/2 13mm,rnd/BOX-36

Product Code :

LBNY-0021-2200080



www.laboratoryinstrumentindia.com

Description :

Type of thread: Non absorbable, synthetic monofilament (Nylon, polyamide).

Length: 75cm.

Initial sterilisation method:

Ethylene oxide gas.

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC 1 (5/0).

Type of needle:

Point: round.

Shape: curved needle, 1/2 circle.

Length: 13mm.

Sterile.

Disposable.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg

Laboratory Instrument India

