



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sut,nonabs,DEC2,need 3/8 18mm,tri/BOX-36

Product Code :

LBNY-0021-220081



www.laboratoryinstrumentindia.com

Description :

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC 2 (3/0).

Type of thread: Non absorbable, synthetic monofilament (Nylon).

Length: Approximately 75cm.

Type of needle:

Point: triangular, conventional cutting

Shape: curved needle, 3/8 circle.

Length: approximately 18mm.

Sterile.

Disposable.

Initial sterilisation method:

Ethylene oxide gas.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg

Laboratory Instrument India

