



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Sut,nonabs,DEC4,need 1/2 30mm,rnd/BOX-36

**Product Code :**

LBNY-0021-220082



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Gauge of thread: Expressed in Decimal gauge DEC number and in the old system of thread numbering: DEC 4 (1).

Type of thread: Non absorbable, synthetic monofilament (Nylon, polyamide).

Length: 75cm.

Type of needle:

Point: round.

Length: 30mm.

Shape: curved needle, 1/2 circle.

Sterile.

Disposable.

Initial sterilisation method:

Ethylene oxide gas.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.030cdm

Estimated weight: 0.010kg

---

# Laboratory Instrument India

