



Email : sales@laboratoryinstrumentindia.com

Product Name :

Clamp,umbilical cord,ster,s.u.

Product Code :

LBNY-0021-2200084



www.laboratoryinstrumentindia.com

Description :

Disposable sterile umbilical cord clamp, safe, effective and easy to apply.

Material: Plastic (polymer).

Grooves all along the length to prevent slip of the umbilical cord and to retain it in the same position.

Finger grip to ensure safe and easy handling.

Safe security lock, with a click to indicate correct locking, to protect against accidental re-opening after clamping.

Sterile.

Single-use.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.16cdm

Estimated weight: 0.004kg



Laboratory Instrument India