



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blanket,survival,220x140cm

Product Code :
LBNY-0021-2200085



www.laboratoryinstrumentindia.com

Description :

Material: Non-stretch polyester film
Lightweight flexible wrap-around blanket
Easy to unfold, strong, tear resistant
Single use, non-sterile
Wind and waterproof, and rot proof
Stays flexible in freezing temperatures
Rectangular size, approximately: 220 x 140cm
One side silver coated, reflects up to 90% of radiated heat
Dimensions:
Length, approximately: 220cm
Width, approximately: 140cm
Thickness, approximately: 12µm

Technical Specification :

Weight/Volume/Dimensions:
Estimated volume: 5cdm
Estimated weight: 0.075kg

Laboratory Instrument India

