



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Blanket,survival,220x140cm

**Product Code :**  
LBNY-0021-2200085



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Non-stretch polyester film  
Lightweight flexible wrap-around blanket  
Easy to unfold, strong, tear resistant  
Single use, non-sterile  
Wind and waterproof, and rot proof  
Stays flexible in freezing temperatures  
Rectangular size, approximately: 220 x 140cm  
One side silver coated, reflects up to 90% of radiated heat  
Dimensions:  
Length, approximately: 220cm  
Width, approximately: 140cm  
Thickness, approximately: 12µm

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 5cdm  
Estimated weight: 0.075kg

---

# Laboratory Instrument India

