



Email : sales@laboratoryinstrumentindia.com

Product Name :

Needle,disp,22G,ster/BOX-100

Product Code :

LBNY-0021-2200087



www.laboratoryinstrumentindia.com

Description :

Material: Needle: Stainless steel

Measurements: External diameter expressed in Gauge and millimetres.

Base and protective cap: medical grade plastic.

Length expressed in inch and millimetres.

Size selected: 22Gx1 1/4" (0.7x30mm).

Sterile

Disposable.

Technical Specification :

Needle,disp,22G,ster/BOX-100



Laboratory Instrument India