



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Needle,disp,25G,ster/BOX-100

**Product Code :**  
LBNY-0021-220088



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Needle: Stainless steel.  
Measurements: External diameter expressed in Gauge and millimetres.  
Base and protective cap: medical grade plastic.  
Length expressed in inch and millimetres.  
Sterile  
Disposable.  
Size selected: 25Gx5/8" (0.5x16mm).

**Technical Specification :**

Needle, disposable, 25Gx5/8" (0.5x16mm), sterile, box of 100.



**Laboratory Instrument India**