



Email : sales@laboratoryinstrumentindia.com

Product Name :

Needle,disp,23G,ster/BOX-100

Product Code :

LBNY-0021-2200089



www.laboratoryinstrumentindia.com

Description :

Material: Needle: Stainless steel conform to ISO 9626-1:1991.
Measurements: External diameter expressed in Gauge and millimetres.
Base and protective cap: medical grade plastic.
Length expressed in inch and millimetres.
Sterile
Disposable.
Size selected: 23Gx1" (0.6x25mm).

Technical Specification :

Needle, disposable, 23Gx1" (0.6x25mm), sterile, box of 100.



Laboratory Instrument India