



Email : sales@laboratoryinstrumentindia.com

Product Name :
Syringe,disp,2ml,w/ndl,21G/BOX-100

Product Code :
LBNY-0021-2200090



www.laboratoryinstrumentindia.com

Description :

Syringe, two pieces: barrel with Luer nozzle and piston/plunger or three pieces: barrel with luer nozzle, stopper and piston/plunger.
Intended use: inject fluids into, or withdraw fluids from, the body
Blister packed.

Technical Specification :

Measurements: External diameter expressed in Gauge and millimeters.
Needle Material: Stainless steel
Base and protective cap: medical grade plastic.
Length expressed in inch and millimeters.
Size selected: 21Gx1 1/2" (0.8x40mm).
Syringe with needle: bi-packed/ detached/mounted needle
Sterile
Disposable.



Laboratory Instrument India

