



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Syringe,disp,2ml,w/ndl,21G/BOX-100

**Product Code :**

LBNY-0021-2200090



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Syringe, two pieces: barrel with Luer nozzle and piston/plunger or three pieces: barrel with luer nozzle, stopper and piston/plunger.

Intended use: inject fluids into, or withdraw fluids from, the body

Blister packed.

**Technical Specification :**

Measurements: External diameter expressed in Gauge and millimeters.

Needle Material: Stainless steel

Base and protective cap: medical grade plastic.

Length expressed in inch and millimeters.

Size selected: 21Gx1 1/2" (0.8x40mm).

Syringe with needle: bi-packed/ detached/mounted needle

Sterile

Disposable.



**Laboratory Instrument India**

