



Email : sales@laboratoryinstrumentindia.com

Product Name :

Syringe,disp,5ml,w/ndl,21G/BOX-100

Product Code :

LBNY-0021-2200091



www.laboratoryinstrumentindia.com

Description :

Syringe, two pieces: barrel with Luer nozzle and piston/plunger or three pieces: barrel with luer nozzle, stopper and piston/plunger.

Intended use: inject fluids into, or withdraw fluids from, the body.

Blister packed.

Technical Specification :

Syringe capacity: 5ml

Syringe Material: medical grade plastic (Polyethylene (PE)/ polypropylene (PP)/ polystyrene (PS).

Needle Material: Stainless steel

Base and protective cap: medical grade plastic.

Barrel Nozzle: The Luer nozzle shall be situated centrally.

Needle Components: Needle with base, and protective cap.

Measurements: External diameter expressed in Gauge and millimeters.

Length expressed in inch and millimeters.

Size selected: 21Gx1 1/2" (0.8x40mm).

Syringe graduation: Graduated scale on the barrel, easy to read. Scale intervals of 0.2 or 0.5ml. Graduation lines of 1.0ml volume increment.

Syringe Barrel: sufficiently transparent to allow easy measurement of the volume contained in the syringe and detection of air bubbles. Length with a maximum usable capacity of at least 10% more than the nominal capacity.



Laboratory Instrument India