



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Syringe,disp,10ml,w/ndl,21G/BOX-100

**Product Code :**

LBNY-0021-2200092



www.laboratoryinstrumentindia.com

**Description :**

Syringe, two pieces: barrel with Luer nozzle and piston or three pieces barrel with luer nozzle, piston and stopper.

Material: Polyethylene (PEF) or polypropylene (PP).

Capacity: 10ml.

Barrel sufficiently transparent to allow easy measurement of the volume contained in the syringe and detection of air bubbles.

Concentric or eccentric luer nozzle.

Graduated scale on the barrel, with scale interval of 0.5 and 1.0ml increment between graduation lines to be numbered.

Material: Needle: stainless steel; base and protecting cap: plastic.

Needle with base (conical Luer type fitting) and protecting cap.

Disposable.

Sterile.

Sterilisation method: Ethylene oxide.

Size selected: 21G (0.80 x 38-40mm).

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 7.917cdm

Estimated weight: 0.866kg



**Laboratory Instrument India**