



Email : sales@laboratoryinstrumentindia.com

Product Name :

Syringe,disp,10ml,ster/BOX-100

Product Code :

LBNY-0021-2200096



www.laboratoryinstrumentindia.com

Description :

Intended use: inject fluids into, or withdraw fluids from, the body.

Syringe, two pieces: barrel with Luer nozzle and piston/plunger or three pieces: barrel with luer nozzle, stopper and piston/plunger.

Blister packed.

Technical Specification :

Syringe capacity: 10ml

Syringe Material: medical grade plastic (Polyethylene (PE)/ polypropylene (PP)/ polystyrene (PS).

Syringe Barrel: sufficiently transparent to allow easy measurement of the volume contained in the syringe and detection of air bubbles. Length with a maximum usable capacity of at least 10% more than the nominal capacity.

Barrel Nozzle: Luer nozzle is centric or acentric.

Syringe graduation: Graduated scale on the barrel, easy to read. Scale intervals of 0.2 or 0.5 or 1.0ml.

Graduation lines of 1.0 or 5.0 ml volume.

Sterile

Syringe: Disposable.

Laboratory Instrument India

