



Email : sales@laboratoryinstrumentindia.com

Product Name :

Syringe,disp,20ml,ster/BOX-100

Product Code :

LBNY-0021-2200097



www.laboratoryinstrumentindia.com

Description :

Syringe, two pieces: barrel with Luer nozzle and piston/plunger or three pieces: barrel with luer nozzle, stopper and piston/plunger.

Intended use: inject fluids into, or withdraw fluids from, the body

Blister packed.

Technical Specification :

Syringe capacity: 20ml

Syringe Material: medical grade plastic (Polyethylene (PE)/ polypropylene (PP)/ polystyrene (PS).

Syringe graduation: Graduated scale on the barrel, easy to read. Scale intervals of 2.0 ml. Graduation lines of 5.0 ml or 10.0 volume.

Barrel Nozzle: Luer nozzle is centric or acentric.

Syringe Barrel: sufficiently transparent to allow easy measurement of the volume contained in the syringe and detection of air bubbles. Length with a maximum usable capacity of at least 10% more than the nominal capacity.

Sterile

Syringe: Disposable.

Laboratory Instrument India

