



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Welding Rod

**Product Code :**  
LBNY-0004-820043



**Description :**

**Welding Rod**

**Technical Specification :**

Welding rod

Electric welding rod: E6013 @3.21 60Kg

Electric welding rod: E6013 A3.21 60IQ

electric welding rod: E7016 A3.2 / 20Kg

Electric welding rod: E70 rc A3,2I20Kg

Aluminum WeldingRod(TIG): A2.4 / 10I\$

Stainless Steel Welding Rod(TIG); T308 A 2.0 /10Kg

Stainless Steel Weiding Rod(TIG): T408 A 2.4 /10Kg

Mild steel welding rod(for TIG): ER70s, @2.4/30Kg

CO2 Wire Solid: ASQIkd /2Roll

CO2 Wire Flux Cord: @1.2(75kg) / 2 Roll

CO2 stainless wire: 01.2 (12.5Kg)/I Roll

CO2 aluminum wtre: 01.2 (7Kg) / I Roll

---

# Laboratory Instrument India

