



Email : sales@laboratoryinstrumentindia.com

Product Name :
Body bag,adult

Product Code :
LBNY-0021-22000101



www.laboratoryinstrumentindia.com

Description :

Body bag for packing and transporting human corpses before burial services.
This type of body bag is not appropriate for the management of contagious diseases.

Technical Specification :

Body bag for handling and transporting human corpses
Special linear reinforced material
Material: PEVA
Chloride and carbon free
Thickness: 200 microns (um)
Tear proof and puncture resistant
100% leak proof seams, entirely heat-sealed, width of 5 mm
Carry capacity: up to 120 Kg
Integrated transparent label pouch 10 x 15 cm (l x w) for placement of identification tag
Dimensions: Unfolded: 220 x 100 cm (l x w) ($\pm 10\%$)
Zipper: High quality black or white zipper, smooth running, full length, U shaped
Single use, disposable
Buried in soil, the body bag will remain solid over a period of minimum 5 years
Large metal or sturdy cloth loops on the zip-runner

Laboratory Instrument India

