



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Body bag,adult

**Product Code :**  
LBNY-0021-22000101



www.laboratoryinstrumentindia.com

**Description :**

Body bag for packing and transporting human corpses before burial services.  
This type of body bag is not appropriate for the management of contagious diseases.

**Technical Specification :**

Body bag for handling and transporting human corpses  
Special linear reinforced material  
Material: PEVA  
Chloride and carbon free  
Thickness: 200 microns (um)  
Tear proof and puncture resistant  
100% leak proof seams, entirely heat-sealed, width of 5 mm  
Carry capacity: up to 120 Kg  
Integrated transparent label pouch 10 x 15 cm (l x w) for placement of identification tag  
Dimensions: Unfolded: 220 x 100 cm (l x w) ( $\pm 10\%$ )  
Zipper: High quality black or white zipper, smooth running, full length, U shaped  
Single use, disposable  
Buried in soil, the body bag will remain solid over a period of minimum 5 years  
Large metal or sturdy cloth loops on the zip-runner

---

# Laboratory Instrument India

