



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Body bag,child

**Product Code :**  
LBNY-0021-22000102



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Body bag Body bag for handling and transporting human corpses before burial services. This type of body bag is not appropriate for the management of contagious diseases (e.g. Ebola haemorrhagic fever, cholera).

**Technical Specification :**

Body bag for handling and transporting human corpses  
Tear proof and puncture resistant  
Single use, disposable  
Material: PEVA  
Special linear reinforced material  
Chloride and carbon free  
Non-porous and non-degradable  
Buried in soil, the body bag will remain solid over a period of minimum 5 years  
100% leak proof seams, entirely heat-sealed, width of 5 mm  
Thickness: 200 microns (um)  
Color: white or black  
Carry capacity: min. 80 kg  
Integrated transparent label pouch 10 x 15 cm (l x w) for placement of identification tag  
High quality black or white zipper, smooth running  
Full length U or J-shaped  
Dimensions:  
Unfolded: 120 x 80 cm (l x w) ( $\pm 10\%$ )

---

Zipper:

Large metal or sturdy cloth loops on the zip-runner

With 25 cm return allowing to show the face of the deceased person.



**Laboratory Instrument India**