



Email : sales@laboratoryinstrumentindia.com

Product Name :
Body bag,child

Product Code :
LBNY-0021-22000102



www.laboratoryinstrumentindia.com

Description :

Body bag Body bag for handling and transporting human corpses before burial services. This type of body bag is not appropriate for the management of contagious diseases (e.g. Ebola haemorrhagic fever, cholera).

Technical Specification :

Body bag for handling and transporting human corpses
Tear proof and puncture resistant
Single use, disposable
Material: PEVA
Special linear reinforced material
Chloride and carbon free
Non-porous and non-degradable
Buried in soil, the body bag will remain solid over a period of minimum 5 years
100% leak proof seams, entirely heat-sealed, width of 5 mm
Thickness: 200 microns (um)
Color: white or black
Carry capacity: min. 80 kg
Integrated transparent label pouch 10 x 15 cm (l x w) for placement of identification tag
High quality black or white zipper, smooth running
Full length U or J-shaped
Dimensions:
Unfolded: 120 x 80 cm (l x w) ($\pm 10\%$)

Zipper:

Large metal or sturdy cloth loops on the zip-runner

With 25 cm return allowing to show the face of the deceased person.



Laboratory Instrument India