



Email : sales@laboratoryinstrumentindia.com

Product Name :

Body bag,infection control,adult

Product Code :

LBNY-0021-22000103



www.laboratoryinstrumentindia.com

Description :

Body bag including absorbent pad used for packing and transporting contagious human corpses before burial services. To be used for management of highly contagious diseases (e.g. Ebola haemorrhagic fever, cholera). Body bag, padded, infection control, adult size.

Technical Specification :

Body bag with absorbent pad prefixed inside of the body bag

Material: PEVA

Chloride and carbon free

100% leak proof seams, entirely heat-sealed, width minimum of 5 mm

Thickness: 300 to 400 microns (um)

Body bag dimensions (with prefixed absorbent pad inside):

Unfolded: 230 x 100 cm (l x w) (\pm 10%)

Carry capacity: min. 120 kg

Special linear reinforced material, minimum: 2 ply

Tear proof and puncture resistant

Non-porous and non-degradable

Color: white

Buried in soil, the body bag will remain solid over a period of minimum 5 years

Minimum of 6 padded carry handles. Non splintering.

Reinforced and integrated into the bag material.

Integrated transparent label pouch 10 x 15 cm (l x w) for placement of identification tag

Body bag is single use, disposable.



Laboratory Instrument India