



Email : sales@laboratoryinstrumentindia.com

Product Name :

Body bag,infection control,child

Product Code :

LBNY-0021-22000104



www.laboratoryinstrumentindia.com

Description :

Body bag including absorbent pad used for packing and transporting contagious human corpses before burial services. To be used for management of highly contagious diseases (e.g. Ebola haemorrhagic fever, cholera).
Body bag, padded, infection control, child size

Technical Specification :

Body bag with absorbent pad prefixed inside of the body bag
Material: PEVA
Special linear reinforced material, minimum: 2 ply
Tear proof and puncture resistant
Non-porous and non-degradable
Chloride and carbon free
Body bag is single use, disposable
Minimum of 6 padded carry handles. Non splintering.
100% leak proof seams, entirely heat-sealed, width minimum of 5 mm
Thickness: 300 to 400 microns (um)
Body bag dimensions (with prefixed absorbent pad inside):
Carry capacity: min. 90 kg
Color: white
Unfolded: min. 150 x 80 cm (l x w) (\pm 10%)
Buried in soil, the body bag will remain solid over a period of minimum 5 years
Integrated transparent label pouch 10 x 15 cm (l x w) for placement of identification tag
Reinforced and integrated into the bag material.



Laboratory Instrument India