



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Coverall, protection, CatIII, Type3b, M

**Product Code :**  
LBNY-0021-22000107



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Disposable liquid-tight biohazard protective coverall, for use in EVD patient isolation units for infection prevention and control against viral penetration.  
After use, the coverall is disposed in a biohazard waste container, collected and destroyed. This is also applicable if the coverall is damaged (perforation, etc.).

**Technical Specification :**

Elasticated hood around face.  
Material: Light weight, do not contain rubber/ latex.  
Protective seams providing barrier equal to fabric.  
Zipper with re-sealable flap protecting leakage through seams.  
Elasticated cuffs and ankles.  
Sleeves with elasticated thumb loop.  
Fabric is Infective agent tested against viral penetration at minimum 1.75kPa  
Non-sterile  
Antistatic treated on both sides.  
Single Use, disposable  
Size: M  
Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.  
Color: White/ yellow/orange.



**Laboratory Instrument India**