



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coverall, protection, CatIII, Type3b, XL

Product Code :
LBNY-0021-22000108



www.laboratoryinstrumentindia.com

Description :

Disposable liquid-tight biohazard protective coverall, for use in EVD patient isolation units for infection prevention and control against viral penetration. After use, the coverall must be disposed in a biohazard waste container, collected and destroyed. This is also applicable if the coverall is damaged.

Technical Specification :

Material: Lightweight, do not contain rubber/ latex.
Protective seams providing barrier equal to fabric.
Zipper with re-sealable flap protecting leakage through seams.
Elasticated hood around face.
Elasticated cuffs and ankles.
Sleeves with elasticated thumb loop.
Fabric is Infective agent tested against viral penetration at minimum 1.75kPa
Antistatic treated on both sides.
Non-sterile
Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.
Color: White/ yellow/orange
Single Use, disposable
Size: XL.



Laboratory Instrument India