



Email : sales@laboratoryinstrumentindia.com

Product Name :

Coverall, protection, CatIII, Type 3b, L

Product Code :

LBNY-0021-22000109



www.laboratoryinstrumentindia.com

Description :

Disposable liquid-tight biohazard protective coverall, for use in EVD patient isolation units for infection prevention and control against viral penetration. After use, the coverall must be disposed in a biohazard waste container, collected and destroyed. This is also applicable if the coverall is damaged.

Technical Specification :

Material: Lightweight, do not contain rubber/ latex.

Antistatic treated on both sides.

Elasticated hood around face.

Elasticated cuffs and ankles.

Sleeves with elasticated thumb loop.

Color: White/ yellow/orange

Fabric is Infective agent tested against viral penetration at minimum 1.75kPa

Non-sterile

Protective seams providing barrier equal to fabric.

Zipper with re-sealable flap protecting leakage through seams.

Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.

Single Use, disposable

Size: L.



Laboratory Instrument India