



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Coverall, protection, CatIII, Type 3b, L

**Product Code :**

LBNY-0021-22000109



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Disposable liquid-tight biohazard protective coverall, for use in EVD patient isolation units for infection prevention and control against viral penetration. After use, the coverall must be disposed in a biohazard waste container, collected and destroyed. This is also applicable if the coverall is damaged.

**Technical Specification :**

Material: Lightweight, do not contain rubber/ latex.

Antistatic treated on both sides.

Elasticated hood around face.

Elasticated cuffs and ankles.

Sleeves with elasticated thumb loop.

Color: White/ yellow/orange

Fabric is Infective agent tested against viral penetration at minimum 1.75kPa

Non-sterile

Protective seams providing barrier equal to fabric.

Zipper with re-sealable flap protecting leakage through seams.

Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.

Single Use, disposable

Size: L.



**Laboratory Instrument India**