



Email : sales@laboratoryinstrumentindia.com

Product Name :
Welding Magnet

Product Code :
LBNY-0004-820044



Description :

Welding Magnet

Technical Specification :

Welding magnet

Safe & Secure: These welding magnets and clamps grip metal and hold it at the angle you need to work safely and efficiently. The welding holders are compact with durable high gloss enamel end plates

Three Angles: Heavy-duty arrow design steel welding magnet square supports material at 45°, 90°, and 135° for easy access. Can be used as double sided holder for welding, assembly and marking work

Multipurpose: The magnetic welding holder clamps connects parts, freeing up your hands for efficient work and improve work safety, and it can also be used as a floater to separate steel plates

Efficiency: Welding Magnet Set Magnetic Tool Holder s ideal welding tools and accessries for being used as holder and positioner in welding, soldering, marking off, pipe installation, making your welding job more efficiently and precisely.

25lbs 30Ea

55lbs 30Ea

75lbs 30Ea

1101bs 30Ea

170lbs 30Ea

Laboratory Instrument India

