



Email : sales@laboratoryinstrumentindia.com

Product Name :

Coverall, protection, CatIII, type 6b, L

Product Code :

LBNY-0021-22000112



www.laboratoryinstrumentindia.com

Description :

Spray/aerosol-penetration resistant, biohazard-protective coverall, for use in isolation units for infection prevention and control against viral penetration.

Coverall, protection, Category III, type 6b, size L.

Technical Specification :

Elasticated hood around face.

Elasticated cuffs and ankles.

Material: Light weight, do not contain rubber/ latex.

Antistatic treated on both sides.

Sleeves with elasticated thumb loop.

Stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.

Stitched or Bound seams.

Zipper with re-sealable flap protecting leakage through seams.

Fabric is Infective agent tested against viral penetration at minimum 1.75kPa (class 2) (or equivalent international standard)

Color: White

Single Use, disposable

Non-sterile

Size: L.



Laboratory Instrument India