



Email : sales@laboratoryinstrumentindia.com

Product Name :

Coverall, protection, Cat III, type 4b, L

Product Code :

LBNY-0021-22000113



www.laboratoryinstrumentindia.com

Description :

Spray-tight biohazard-protective coverall, for use in EVD patient-isolation units for infection prevention and control against viral penetration.

Technical Specification :

Elasticated hood around face.

Elasticated cuffs and ankles.

Material: Lightweight, do not contain rubber/ latex.

Antistatic treated on both sides. Fabric is Infective agent tested against viral penetration at minimum 1.75kPa (Class 2) (or equivalent international standard)

Non-sterile

Sleeves with elasticated thumb loop.

Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.

Color: White

Protective seams providing barrier equal to fabric.

Zipper with re-sealable flap protecting leakage through seams.

Single Use, disposable

Size : L.

Laboratory Instrument India

