



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coverall, protection, Cat III, type 4b, XL

Product Code :
LBNY-0021-22000114



www.laboratoryinstrumentindia.com

Description :

Spray-tight biohazard-protective coverall, for use in EVD patient-isolation units for infection prevention and control against viral penetration.

Technical Specification :

Protective seams providing barrier equal to fabric.
Zipper with re-sealable flap protecting leakage through seams.
Elasticated hood around face.
Sleeves with elasticated thumb loop.
Elasticated cuffs and ankles.
Material: Lightweight, do not contain rubber/ latex.
Antistatic treated on both sides. Fabric is Infective agent tested against viral penetration at minimum 1.75kPa (Class 2) (or equivalent international standard)
Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.
Color: White
Single Use, disposable
Non-sterile
Size : XL.

Laboratory Instrument India

