



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Coverall, protection, CatIII, type 4b, M

**Product Code :**  
LBNY-0021-22000115



www.laboratoryinstrumentindia.com

**Description :**

Spray-tight biohazard-protective coverall, for use in EVD patient-isolation units for infection prevention and control against viral penetration.

**Technical Specification :**

Elasticated hood around face.  
Elasticated cuffs and ankles.  
Material: Lightweight, do not contain rubber/ latex.  
Sleeves with elasticated thumb loop.  
Color: White  
Antistatic treated on both sides. Fabric is Infective agent tested against viral penetration at minimum 1.75kPa (Class 2) (or equivalent international standard)  
Protective seams providing barrier equal to fabric.  
Zipper with re-sealable flap protecting leakage through seams.  
Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.  
Non-sterile  
Single Use, disposable  
Size : M

---

# Laboratory Instrument India

