



Email : sales@laboratoryinstrumentindia.com

Product Name :

Apron,protection,plastic,reusable

Product Code :

LBNY-0021-22000117



www.laboratoryinstrumentindia.com

Description :

Reusable straight sleeveless protective apron, for use in healthcare settings.

Material: Opaque or translucent high quality plastic material.

Water and body fluids resistant, can be disinfected with chlorine solution (0.05% concentration).

Back fastening band and neckband should be strong and not detachable.

Moisture-proof and stain resistant.

Straight apron with bib, back fastening and neckband.

Reusable

Non-sterile

Size selected: Adult size.

Length: 120cm (+/-5cm) from top of the bib to lower edge of the apron, total length: 145cm (+/-5cm).

Width: 80cm (+/-10cm).

Thickness: 150-300 microns.

Technical Specification :

Reusable straight sleeveless protective apron, for use in healthcare settings.



Laboratory Instrument India

