



Email : sales@laboratoryinstrumentindia.com

Product Name :
Non-pneumatic Anti-shock Garment,sizeM/L

Product Code :
LBNY-0021-22000118



www.laboratoryinstrumentindia.com

Description :

A medium/large sized, non-pneumatic anti-shock garment (NASG) to be used for managing uncontrollable postpartum haemorrhage (PPH) and to keep women alive until they can get the treatment they need.

Technical Specification :

Material: Made of a lightweight neoprene garment and Velcro.
Lightweight, flexible and comfortable for the wearer.
Reusable
No metal parts: safe for X-rays and MRIs.



Laboratory Instrument India