



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Bootcover,antiskid,elasticated,pair

**Product Code :**

LBNY-0021-22000120



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Overboot cover to provide partial body protection to prevent cross-contamination between the boot and the environment.

Single use only.

Heavy duty over boot cover, puncture and abrasion resistant. Single use, pair.

**Technical Specification :**

Material: polypropylene, coated, non-woven

Seam type: bound seams

Antistatic treated

Anti-slip sole

Heavy duty over boot cover, puncture- and abrasion resistant

Total length: sole till upper-calf level: 45 cm

Sole thickness: 0.55 mm

Size: one size (fits EU sizes 42-46)

Fits either foot (ambidextrous)

Liquid impermeable and repellent (protection against bio-hazardous liquids)

Elasticated opening

With fixation ties

Single use, disposable

Colour: white

Unit is pair.



**Laboratory Instrument India**