



Email : sales@laboratoryinstrumentindia.com

Product Name :

Apron,protection,plastic,disp

Product Code :

LBNY-0021-22000121



www.laboratoryinstrumentindia.com

Description :

To be used in healthcare facilities by personnel performing medical/obstetrical/surgical procedures or cleaning with high risk of contamination by body fluids projection.

Technical Specification :

Straight, sleeveless protective apron, for use in healthcare settings
Seamless
Material: polyethylene (PE) or biodegradable or compostable material
Liquid proof and stain resistant
Both back- and neck-band can be adjusted/fastened
Comfortable to wear, apron has back- and neck-band strips attached (4 in total)
Single-use
Can resist water and disinfectant (ethanol 70% and chlorine solution 0.5%)
Size: 85 x 145 cm (w x l) (+/- 15%)
Thickness: min: 50 um
Color: white.



Laboratory Instrument India

