



Email : sales@laboratoryinstrumentindia.com

Product Name :

Cap,surgical,bouffant,non-woven,BOX-100

Product Code :

LBNY-0021-22000125



www.laboratoryinstrumentindia.com

Description :

To be worn by surgical staff during surgical procedures or in other health settings to protect both patient and medical staff from cross-contamination. An elasticated cap to completely cover the wearers hair.

Technical Specification :

Material: nonwoven polypropylene
Nonwoven surgical cap
Round, bouffant surgical cap, elasticated
Non-sterile
Single use
One size: Adult, standard size
Colors: blue or green.



Laboratory Instrument India